

Depyrogenation Study on Primary Packaging

Michael Kracklauer¹, Werner Iländer², Peter Alexander Kitschmann², Maike Piehler¹, Johannes Reich¹



1 Microcoat Biotechnologie GmbH

Am Neuland 3, 82347 Bernried am Starnberger See, Germany
Tel.: +49 81 58 99 81 0
Email: info@microcoat.de
www.microcoat.de



2 Bausch + Ströbel Maschinenfabrik

Ilshofen GmbH + Co. KG

Parkstraße 1, 74532 Ilshofen, Germany
Tel.: +49 7904 7010
Email: info@bausch-stroebel.de
www.bausch-stroebel.com

Abstract


Endotoxin (LPS), if present in the bloodstream, is already in small doses a severe health risk. Thus depyrogenation of primary packing is important to ensure patient safety. Dry-heat depyrogenation is the primary method used for endotoxin inactivation. A depyrogenation process should demonstrate at least 99.9 % (3-log) endotoxin reduction. Common processes used in the pharmaceutical industry to prepare vials before filling with parenteralia often include sonication, washing and heating. Sonication and washing are reducing the number of particles in the glass vials, but are these vial treatments able to reduce the contamination with endotoxin as well or is heating absolutely necessary? Very few data is available about endotoxin reduction in vials at different glass depyrogenation treatments.

In this study a systematic approach was used to test 6 different combinations of vial treatments and their influence on endotoxin contamination of 5 different glass vials (n = 9-10).



Method and Materials

Workflow




→ Purified LPS from *E.coli* O113

- Drying: 1000 EU/vial, over night
- Vial Treatments: 6 different methods (0, a-e), negative controls (NC) without addition of endotoxin and undried LPS (R) was used as controls (each n = 9-10)
- Reconstitution: in 0.2 % SDS for 2.5 h including vortexing
- Sample dilution in LRW: 1:100, 1:1000
- Detection: rFC Method (Endozyme II)

5 different glass vials used for the study

- Endograde Glass Test Tube (Hyglos)
- ISO 2R Clear Glass Vial (Nuova Ompi)
- ISO 2R Flacon (Gerresheimer Chalon)
- ISO 2R Fiolar clear StandardLine (Schott)
- ISO 2R Fiolar clear TopLine (Schott)

Overview – 6 different vial treatments



Sonicator

1 min 70 °C
240 W

Cleaning machine


3 x wasching
water temperature 80 °C

Heat tunnel

15 min
320 °C

Results

Endotoxin concentration after vial treatment Mean ± SD (n = 9-10)



Log reduction of endotoxin per vial

Vial treatment	Endograde	Ompi	Gerresheimer	Schott (SL)	Schott (TL)
a	0.93	0.55	0.88	0.57	1.19
b	1.18	0.99	1.05	0.72	1.32
c	> 3.39	> 3.41	> 3.30	> 3.9	> 3.49
d	> 3.39	> 3.41	> 3.30	> 3.29	> 3.49
e	> 3.39	> 3.41	> 3.30	> 3.29	> 3.49
NC	> 3.39	> 3.41	> 3.30	> 3.29	> 3.49

- LPS concentration is not influenced by drying method (0 vs. R)
- 0.57-log to 1.32-log LPS reduction for treatment a and b
- > 3-log reduction for treatments c, d and e (equal NC)

Conclusion

- Endotoxin in different glass vials behaves comparable with and without treatments
- For reconstitution of endotoxin a dedicated sample preparation (0.2 % SDS) was applied
- Washing and sonicating is not sufficient to reduce endotoxin (reduction ≤ 1.32-log)

→ **A 3-log endotoxin reduction was only possible after heat treatment**

References

- Tech Report #3 Validation of Dry Heat Processes Used for Sterilization and Depyrogenation, PDA
- USP <1211> Dry-Heat Sterilization/Depyrogenation

